

Microbiome Testing Line

Clinical Product Portfolio

Comprehensive real-time PCR-based molecular diagnostics for quantitative microbiota assessment across gynecological, pediatric, and urological applications



FEMOBIOME II LINE

Women's Urogenital Microbiota

- 3 test variants: FEMOBIOME II (16 tubes), PrimaScreen (8 tubes), SecundaScreen (8 tubes)
- Normobiota, anaerobes, aerobes, mycoplasmas, STIs, HPV, herpesviruses
- Vaginal and cervical sampling incl. self-collection
- Visual traffic light reporting system

ENTEROBIOME KIDS

Children's Gut Microbiota

- Covers 99.9% of prokaryotic mass in children's feces
- 6 bacterial phyla + Euryarchaeota + Candida
- Pathogenicity markers: mecA, srr2, tcdA, tcdB
- Normocenosis/dysbiosis status with diversity estimation

ANDROBIOME

Male Urogenital Microbiota

- 2 variants: Androbiome (full) and Androbiome Screen
- No analogues worldwide
- Normal, opportunistic microbiota and STI detection
- 6 biomaterial types supported

Expanded line of tests for comprehensive women's urogenital microbiota profiling with real-time PCR. Features optimized pathogen panels, quantitative GBS assessment for pregnancy, and viral detection including HPV and herpesviruses.

What's New in FEMOBIOME II

- Expanded and optimized microorganism detection panel
- More informative normobiota evaluation across age groups
- Quantitative Group B Streptococcus assessment for pregnancy planning
- Obligate pathogen control for asymptomatic infections
- HPV detection: types 16, 18, 45 + 11 additional high-risk types
- Herpesvirus detection: HSV-1, HSV-2, CMV

Three Test Variants

FEMOBIOME II

Universal test — 16 tubes. Most comprehensive capabilities for full microbiota profiling.

PrimaScreen

Shortened version — 8 tubes. Focused on bacterial vaginosis diagnostics.

SecundaScreen

Complex test — 8 tubes. Infectious and inflammatory urogenital disease diagnostics.

FEMOBIOME II — Detection Panel Comparison

Comprehensive comparison of microorganisms detected by each test variant

Category	Organism / Marker	FB II	Prima	Secunda
Control	Human DNA / Total bacterial load	✓	✓	✓
Normobiota	Lactobacillus spp., L. crispatus, L. jensenii, L. gasseri, L. iners, Bifidobacterium	✓	Partial	Partial
Aerobes	Staphylococcus, Streptococcus, S. agalactiae, Enterobacteriaceae, Enterococcus, Haemophilus	✓	—	✓
Anaerobes	Gardnerella, Fannyhessea, Mobiluncus, Anaerococcus, Peptostreptococcus, Bacteroides/Prevotella, Sneathia, Megasphaera, BVAB1-3	✓	Key only	✓
Mycoplasmas	U. urealyticum, U. parvum, M. hominis	✓	✓	✓
Yeast	Candida spp., Candida albicans	✓	✓	✓
STIs	C. trachomatis, M. genitalium, N. gonorrhoeae, T. vaginalis	✓	—	✓
Herpes	HSV-1, HSV-2, CMV	✓	—	✓
HPV	HPV 16, 18, 45 + 11 high-risk types	✓	—	—

Clinical Indications

Infectious Diseases

Inflammatory symptoms → FEMOBIOME II; BV symptoms → FEMOBIOME II or SecundaScreen

Asymptomatic Cases

Preparation for pelvic surgery; Preventive screening → FEMOBIOME II or PrimaScreen

Treatment Monitoring

Post-treatment assessment of infectious diseases and BV

Obstetrics

Pregnancy planning and ART; Prophylactic exams during pregnancy; Infertility and miscarriage workup → FEMOBIOME II

Biomaterial & Preparation

Sample Types

- Vaginal swab (specialist or self-collected)
- Cervical swab in transport medium
- Compatible with liquid cytology media (PreservCyt, BD SurePath, EASYPREP, CellIPrep)

Patient Preparation

- 24h: No syringing, tampons, antiseptics, transvaginal ultrasound
- 48h: No colposcopy
- 72h: No unprotected sexual contacts
- 2 weeks after antibiotic/probiotic therapy
- 4 weeks after STI therapy

Traffic Light Indicator System

State	Pathogens Detected	No Pathogens Viruses/Myco/Candida $\geq 10^4$	No Pathogens No Significant Findings
Eubiosis	RED	YELLOW	GREEN
Moderate Dysbiosis	RED	YELLOW	YELLOW
Severe Dysbiosis	RED	RED	RED

Mycoplasmas: *U. urealyticum*, *U. parvum*, *M. hominis*. Candida: *C. spp.*, *C. albicans*. Thresholds: $\geq 10^4$ GE/ml.

Report Components

- Indicator icon — visual summary combining microbiota state and pathogen presence
- Linear histogram — proportions of normobiota, anaerobes, aerobes, and mycoplasmas with color-coded quantity levels
- Columnar histogram — detailed distribution of conditionally pathogenic microorganisms within each group

Result Categories

Personalized comments and infographics for each patient

-  Eubiosis — no pathogens, no significant viruses/mycoplasmas/Candida
-  Warning — moderate dysbiosis or elevated viruses/mycoplasmas/Candida
-  Alarm — severe dysbiosis or obligate pathogen detected

✓ Eubiosis

Normal microbiota dominance

Conclusion

Microbiota status — eubiosis: normobiota dominance, *Lactobacillus* spp. relative proportion 99.8% (*L. crispatus*, *L. iners*, *L. jensenii*).

No pathogens detected.

Biomaterial Control Parameters

Human genomic DNA: 4.9 Ig GE/ml (Ref: ≥ 3.5)

Total bacterial load: 6.2 Ig GE/ml (Ref: ≥ 4.0)

Pathogens — All Not Detected

SEXUALLY TRANSMITTED INFECTIONS

Chlamydia trachomatis — not detected
Mycoplasma genitalium — not detected
Neisseria gonorrhoeae — not detected
Trichomonas vaginalis — not detected

HERPESVIRUSES

HSV-1 — not detected
HSV-2 — not detected
CMV — not detected

HUMAN PAPILLOMAVIRUSES (HPV)

HPV 16 — not detected
HPV 18 — not detected
HPV 45 — not detected
HPV 31/33/35/39/51/52/56/58/59/66/68 — not detected

Microbiota Composition

Example: Eubiosis case — vaginal swab

Normobiota 99.8%

Threshold: $\geq 80\%$ normobiota = eubiosis

Micro organism	Ig GE/ml	%
NORMOBIOTA		
<i>Lactobacillus</i> spp.	6.2	99.8%
<i>Bifidobacterium</i> spp.	1.7	<1%
AEROBES		
<i>Staphylococcus</i> spp.	—	—
<i>Streptococcus</i> spp.	—	—
Enterobacteriaceae	—	—
<i>Enterococcus</i> spp.	2.4	<1%
<i>Haemophilus</i> spp.	—	—
ANAEROBES		
<i>Gardnerella vaginalis</i>	—	—
<i>Fannyhessea vaginae</i>	—	—
<i>Mobiluncus</i> spp.	—	—
<i>Anaerococcus</i> spp.	2.9	<1%
<i>Peptostreptococcus</i> spp.	2.6	<1%
<i>Bacteroides</i> / <i>Prevotella</i>	3.3	<1%
<i>Sneathia</i> / <i>Leptotrichia</i> / <i>Fusobacterium</i>	—	—
<i>Megasphaera</i> / <i>Veillonella</i> / <i>Dialister</i>	2.4	<1%
BVAB1 / BVAB2 / BVAB3	—	—
MYCOPLASMAS		
<i>Ureaplasma urealyticum</i>	—	—
<i>Ureaplasma parvum</i>	<4.0	<1%
<i>Mycoplasma hominis</i>	—	—
YEAST		
<i>Candida</i> spp.	—	—
<i>Candida albicans</i>	—	—

Kit Specifications

- Analysis time: from 2.5 hours (with sample pretreatment)
- PCR formats: Automated (384-well) and manual (96-well)
- Package S: 12 tests (24 for PrimaScreen/SecundaScreen)
- Package A / A-TL: 24 tests (48 for PrimaScreen/SecundaScreen)
- Storage: +2 to +8°C (packages S, A, A-TL)
- Storage: -18 to -22°C (TechnoTaq MAX polymerase)
- Shelf life: 12 months
- Transport medium: STOR-F, STOR-M

Compatible Platforms

Extraction Platforms

- Allsheng Auto-Pure 96
- KingFisher Flex (PREP-MB-RAPID II)

Liquid Cytology Media

- PreservCyt (Hologic Inc., USA)
- BD SurePath (Becton, Dickinson, USA)
- EASYPREP (YD Diagnostics, South Korea)
- CellPrep (CP Biodyne, South Korea)
- Cell Preservative Solution (Human Lituo, China)

Sample Pretreatment

- PREP-PK for scrapes in BD SurePath and CellPrep media

Reagent kit for quantitative analysis of gut microbiota composition in children's stool by real-time PCR. Detects microbial taxa covering 99.9% of total prokaryotic mass.

Detection Scope

- Actinomycetota: Bifidobacterium spp. (child and adult species), Coriobacteria
- Bacillota: Cl. leptum gr., F. prausnitzii, Lachnospiraceae, Lactobacillaceae, Streptococcus spp.
- Bacteroidota: Bacteroides, Prevotella, Parabacteroides, Alistipes, Butyrivimonas spp.
- Verrucomicrobiota: Akkermansia muciniphila
- Pseudomonadota: Enterobacterales
- Euryarchaeota: Methanobacteriaceae
- Candida spp. (yeast fungi)

Key Capabilities

- Normocenosis/dysbiosis status determination
- Normal microbiota diversity estimation
- Bifidobacteria species-level identification (B. bifidum, B. longum, B. infantis, B. breve, B. adolescentis)
- Bacillota/Bacteroidota ratio calculation
- Pathogenicity and resistance markers: mecA (methicillin resistance), srr2 (S. agalactiae invasiveness), tcdA/tcdB (C. difficile toxins)

Kit Features

- Multiplex format with internal control
- Automatic results generation via RealTime_PCR software
- Preset templates for major PCR instruments
- Biomaterial: feces including meconium

Four Phases of Gut Microbiota Formation in Children

Phase 1: Birth

Streptococcus spp., Enterobacterales

Phase 2: 14 Days

Bifidobacterium, Lactobacillaceae, Clostridium,
Bacteroidota

Phase 3: Complementary Foods

Bacteroidota, Bacillota, Verrucomicrobiota,
Actinomycetota

Phase 4: Weaning

Stabilization — adult-like composition

Factors Affecting Microbiota Formation

- **Type of birth:** Natural birth → stable composition; C-section → delayed colonization, more *C. difficile*, *Enterococcus*, *Klebsiella*, *Veillonella*
- **Feeding type:** Breastfed → *Bifidobacterium* dominant (*B. breve*, *B. bifidum*, *B. longum*); Formula-fed → *Bacteroides* and *Clostridium* dominant
- **Gestational age:** Preterm → delayed anaerobic colonization, higher Enterobacteriaceae and opportunistic organisms
- **Antibiotic use:** Reduces diversity, increases drug-resistant Enterobacteriaceae and *Clostridium*, decreases Bifidobacteriaceae

A species-rich gut ecosystem is more resilient to adverse stresses. Reduced diversity is associated with disease presence.

Normal Microbiota

Actinomycetota

- Bifidobacterium spp. (child species: B. bifidum, B. longum, B. infantis, B. breve)
- Bifidobacterium spp. (adult species: B. adolescentis, B. catenulatum, B. animalis)
- Coriobacteria

Bacillota

- Cl. leptum gr. (cluster IV)
- F. prausnitzii
- Lachnospiraceae (cluster XIVa)
- Lactobacillaceae
- Lactococcus lactis
- Streptococcus spp.
- Veillonellaceae
- Erysipelotrichaceae

Bacteroidota

- Bacteroides spp.
- Prevotella spp.
- Parabacteroides spp.
- Alistipes spp.
- Butyrivibrio spp.

Other

- Akkermansia muciniphila (Verrucomicrobiota)
- Blautia spp. (acetogens)
- Desulfovibrio spp. (sulfate-reducing)
- Methanobacteriaceae (methanogens)

Opportunistic & Markers

Pathobionts

- C. difficile (tcdA, tcdB markers)
- C. perfringens gr.
- Enterobacterales
- Enterococcus spp.
- S. aureus (mecA marker)
- S. agalactiae (srr2 marker)

Fungi

- Candida spp.
- Candida albicans

Key Clinical Markers

- **Bacillota/Bacteroidota ratio:** Shift toward Bacteroidota may indicate inflammatory process
- **F. prausnitzii:** Anti-inflammatory biomarker. Promotes IL-10/IL-12 secretion, inhibits IL-8. Decreased in inflammatory bowel disease
- **Akkermansia muciniphila:** Main mucin metabolizer. Enhances barrier function. Produces propionate and acetate for butyrate production
- **Bifidobacterium species:** Child species (B. bifidum, B. longum, B. infantis, B. breve) ferment breast milk oligosaccharides, activate T-regulatory cells. Adult species ferment plant polysaccharides

Disease Associations & Markers

Pathogenicity Markers

- **mecA** — methicillin resistance in Staphylococcus; strains resistant to methicillin
- **srr2** — S. agalactiae invasiveness; promotes bloodstream and CNS invasion, causing meningitis
- **tcdA / tcdB** — C. difficile enterotoxins A and B; cause acute inflammation, intestinal barrier damage, fluid influx

Disease Associations

- IBD (UC, Crohn's): ↑ Bacteroides/Prevotella, ↓ F. prausnitzii
- Celiac disease: ↓ Lactobacillus/Bifidobacteria, ↑ pro-inflammatory Veillonellaceae
- Necrotizing enterocolitis (preterm): associated with C. perfringens and Enterobacterales
- C. difficile carriage: up to 61% in children <1 year; normally eliminated by ~2 years

Result Form Components

1. Normocenosis / Dysbiosis Status

Determined by relative abundance of normobiota vs. opportunistic microbiota in TBM. Classification guides clinical decision-making.

Normocenosis

Moderate Dysbiosis

Severe Dysbiosis

2. Diversity Estimation

Normal microbiota diversity across 6 phyla: Actinomycetota, Bacillota, Bacteroidota, Pseudomonadota, Fusobacteriota, Verrucomicrobiota. Lower diversity correlates with disease risk.

3. Bacillota / Bacteroidota Ratio

Key gut balance indicator. Shift toward Bacteroidota may indicate inflammatory process. Associated with metabolic disorders, obesity, and IBD.

4. Quantitative Assessment

Bifidobacteria/lactobacilli counts with species-level identification. Opportunistic microbiota and *Candida* quantification. Includes *F. prausnitzii* and *Akkermansia muciniphila* as anti-inflammatory biomarkers.

Clinical Associations

IBD: ↑ Bacteroides/Prevotella, ↓ *F. prausnitzii*. Celiac: ↓ Lactobacillus/Bifidobacteria. NEC: *C. perfringens*, Enterobacterales. *C. difficile* carriage patterns tracked.

Pathogenicity Markers

mecA

Methicillin-resistant *Staphylococcus* detection

srr2

Str. agalactiae invasiveness; promotes meningitis

tcdA

C. difficile enterotoxin A

tcdB

C. difficile enterotoxin B (cytotoxin)

Automated Report Generation

Results form auto-generated with recommended Real-Time PCR instruments and RealTime_PCR software. Preset templates configure settings and calculate results automatically. No manual data entry required.

Report Includes

- Quantitative normobiota estimation by phylum
- Bifidobacteria species identification
- Opportunistic microbiota and *Candida* levels
- Antibiotic resistance markers (mecA)
- Virulence markers (srr2, tcdA, tcdB)
- Bacillota/Bacteroidota ratio calculation
- Overall dysbiosis classification

Comprehensive real-time PCR-based assay for male urogenital tract microbiota assessment and STI detection. No analogues worldwide. Hands-on time from 3 hours.

Androbiome (Full)

Detailed microbiome structure assessment

- Normal flora: Staphylococcus, Streptococcus, Corynebacterium, Lactobacillus spp.
- BV-associated: Gardnerella, Megasphaera, Sneathia, Ureaplasma, Mycoplasma, Atopobium
- Anaerobes: Bacteroides, Anaerococcus, Peptostreptococcus, Eubacterium spp.
- Other OM: Haemophilus, Pseudomonas/Ralstonia/Burkholderia, Enterobacteriaceae/Enterococcus
- Yeast: Candida spp.
- STIs: C. trachomatis, M. genitalium, N. gonorrhoeae, T. vaginalis

Androbiome Screen

Screening and acute diagnostics

- Shortened panel for rapid screening
- Normal: Staphylococcus, Streptococcus, Corynebacterium
- Opportunistic: Gardnerella, Ureaplasma, Mycoplasma, Enterobacteriaceae/Enterococcus, Candida
- Full STI panel: C. trachomatis, M. genitalium, N. gonorrhoeae, T. vaginalis
- Ideal for acute urogenital tract diseases

When to Use Which

- Screen: Screening, differential diagnosis of acute forms
- Full: Chronic diseases, ineffective treatment, reproductive function, IVF preparation

ANDROBIOME — Detection Panel Comparison

Complete comparison of organisms detected by Androbiome and Androbiome Screen

Category	Organism	Full	Screen
Normal	Staphylococcus spp., Streptococcus spp., Corynebacterium spp.	✓	✓
Transit	Lactobacillus spp.	✓	✓
BV-associated	Gardnerella, Megasphaera/Veillonella/Dialister, Sneathia/Leptotrichia/Fusobacterium, Atopobium cluster	✓	—
Mycoplasmas	Ureaplasma urealyticum, U. parvum, Mycoplasma hominis	✓	✓
Anaerobes	Bacteroides/Porphyromonas/Prevotella, Anaerococcus, Peptostreptococcus/Parvimonas, Eubacterium	✓	—
Other OM	Haemophilus spp., Pseudomonas/Ralstonia/Burkholderia, Enterobacteriaceae/Enterococcus	✓	Partial
Yeast	Candida spp.	✓	✓
STIs	C. trachomatis, M. genitalium, N. gonorrhoeae, T. vaginalis	✓	✓

Clinical Indications

- **Androbiome Screen:** Screening and differential diagnosis of acute urogenital diseases (urethritis, balanoposthitis)
- **Androbiome (Full):** Chronic forms, ineffective treatment, reproductive function assessment, IVF preparation
- **Couple examination:** Combined with FEMOBIOME for infertility workup

Preparation Requirements

- Sexual abstinence or barrier contraception for 3 days
- Avoid antiseptics and antibacterial soaps
- Avoid urinating 1.5-2 hours before urethral sampling
- Test can be prescribed during antibiotic therapy for monitoring

Biomaterial Options

- **Epithelial scrapes:** Glans penis, foreskin, urethra (urethritis, balanoposthitis)
- **Ejaculate:** Epididymitis, prostatitis, asymptomatic STIs, male infertility
- **Prostatic fluid:** After prostate massage (acute prostatitis: massage prohibited)
- **Residual urine post-massage:** Alternative to prostatic fluid collection
- **First-portion urine:** STI detection only (acute inflammation with pain)
- **Biopsy samples:** Prostatic tissue in sterile saline

Storage: 2-8°C up to 24 hours. Frozen: -18 to -22°C up to 1 month. Multiple PCR tests can use a single sample.

4-Step Interpretation

- Assess obligate pathogens (STIs) and Candida spp. — reported even with insufficient biomaterial
- Control biomaterial compliance: Human DNA $>10^3$, TBM $>10^4$. Both below threshold → insufficient sample
- Assess balance: normal vs opportunistic microbiota. Normal = 2+ genera of commensals at TBM $>10^5$
- Generate automatic conclusion: pathogen presence, Candida status, microbiome structure

Transit Microbiota (Lactobacillus)

- Marker of female partner's microbiota in male tract
- If $>10\%$ of TBM → microbiome composition NOT analyzed
- Recommend: repeat with barrier contraception or 3-day abstinence

Key Indicators

Control Parameters

- Human DNA: shows amount of human genomic material in sample
- Total Bacterial Mass (TBM): total bacteria for proportion calculations
- Total Microorganisms Detected (TMD): mathematical parameter for group ratios
- Candida and T. vaginalis NOT included in TBM

Result Categories

- Normal composition — commensals (Staphylococcus, Streptococcus, Corynebacterium) dominant
- Moderate disturbance — decreased normal microbiota proportion
- Severe disturbance — opportunistic predominance or STI detected

Couple Examination

Combined ANDROBIOME + FEMOBIOME for comprehensive partner evaluation in infertility workup

Real-Time PCR vs. Alternative Methods

vs. Bacterial Culture

- 50-80% of gut microorganisms are difficult to culture
- No dependency on organism viability
- Detects anaerobes that culture methods miss
- Faster turnaround time
- More comprehensive coverage

vs. NGS Sequencing

- Results are clinically actionable without specialized interpretation
- More cost-effective for routine diagnostics
- Standardized and reproducible
- Suitable as routine clinical method
- No bioinformatics expertise required

Key Platform Advantages

- All-in-one: from sample collection to report interpretation
- Automated (384-well) and manual (96-well) PCR formats
- Fast nucleic acid extraction
- Multiplex format for efficiency
- Quantitative results with automatic interpretation

Extraction & PCR Platforms

DNA Extraction

- Allsheng Auto-Pure 96
- KingFisher Flex (PREP-MB-RAPID II)
- PREP-MB MAX (for ejaculate)
- PREP-NA PLUS, PREP-GS PLUS (alternatives)

PCR Setup

- Automated 384-well format
- Manual 96-well format
- Preset templates for major real-time PCR instruments
- Automatic results generation via RealTime_PCR software

Storage & Logistics

Kit Storage

- +2 to +8°C: packages S, A, A-TL
- 18 to -22°C: TechnoTaq MAX polymerase (package A)
- Shelf life: 12 months for all kits

Sample Storage

- 2-8°C: up to 24 hours
- 18 to -22°C: up to 1 month (frozen)
- Transport media: STOR-F, STOR-M

Compatible Liquid Cytology Media

- PreservCyt (Hologic Inc., USA)
- BD SurePath (Becton, Dickinson, USA)
- EASYPREP (YD Diagnostics, South Korea)
- CellPrep (CP Biodyne, South Korea)



Microbiome Testing Solutions

Comprehensive molecular diagnostics for clinical microbiota assessment

FEMOBIOME II

Women's urogenital microbiota

ENTEROBIOME KIDS

Children's gut microbiota

ANDROBIOME

Male urogenital microbiota

For inquiries and orders:

www.genomer.com.tr

info@genomer.com.tr