

PRODUCT CATALOGUE

Molecular Diagnostics Solutions

2026 EDITION

Genomer Molecular Diagnostics

Genomer Biyoteknoloji Arge San. ve Tic. Ltd. Şti. — delivering advanced real-time PCR-based molecular diagnostic solutions for microbiome assessment, pathogen detection, and HPV screening.



Real-time PCR

Gold-standard multiplex qPCR technology



Quantitative Results

Precise quantification with clinical interpretation



Validated Workflows

End-to-end solutions from sample to report

Table of Contents

Microbiome Portfolio Overview	3
FEMOBIOME® II Line	4–6
ANDROBIOME® Line	7
ENTEROBIOME® Kids	8
HPV Testing Overview	9
Genomer HPV Screen	10
Genomer HPV Quant 15	11
Genomer HPV Quant 21	12

Microbiome Portfolio Overview

Comprehensive real-time PCR-based microbiome assessment kits for clinical diagnostics.

FEMOBIOME® II Line

Line of tests for women microbiota assessment

3 test kits 8–16 tubes per kit Vaginal & cervical swabs

ANDROBIOME® Line

Comprehensive real-time PCR-based assay for the assessment of the state of the male urogenital tract microbiota and detection of obligate pathogens — causative agents of sexually transmitted infections (STIs).

2 test kits Urethral scrapes, ejaculate, urine, prostate fluid

ENTEROBIOME® Kids

Reagent kit for the analysis of gut microbiota composition in stool samples from children by real-time PCR.

1 test kit 12 tests per kit Stool samples

FEMOBIOME® II Line

FEMOBIOME®II – universal test for assessment women microbiota state with the most comprehensive capabilities.

FEMOBIOME®PrimaScreen – a shortened version of FEMOBIOME®II for diagnostics of bacterial vaginosis.

FEMOBIOME®SecundaScreen – a complex test for diagnostics of infectious and inflammatory urogenital diseases.

	FB®II	Prima Screen	Secunda Screen
Number of tubes	16	8	8
Indicators			
Human genomic DNA	●	●	●
Total bacterial load	●	●	●
Normobiota			
Lactobacillus spp.	●	●	●
<i>L. crispatus</i>	●		
<i>L. jensenii</i> / <i>L. mulieris</i>	●		
<i>L. gasseri</i> / <i>L. paragasseri</i>	●		
<i>L. iners</i>	●	●	●
<i>L. non-iners</i>		●	●
Bifidobacterium spp.	●		
Aerobes			
Staphylococcus spp.	●		●
Streptococcus spp.	●		●
<i>Streptococcus agalactiae</i>	●		
Enterobacteriaceae	●		●
Enterococcus spp.	●		●
Haemophilus spp.	●		●
Anaerobes (BV-associated)			
Gardnerella vaginalis	●	●	●
Fannyhessea vaginae (Atopobium vaginae)	●	●	●
Mobiluncus spp.	●	●	●
Anaerococcus spp.	●		●
Peptostreptococcus spp.	●		●
Bacteroides spp. / Porphyromonas spp. / Prevotella spp.	●		●
Sneathia spp. / Leptotrichia spp. / Fusobacterium spp.	●		●
Megasphaera spp. / Veillonella spp. / Dialister spp.	●		●
BVAB1 / BVAB2 / BVAB3	●		●

	FB® II	Prima Screen	Secunda Screen
Mycoplasmas			
Ureaplasma urealyticum	●	●	●
Ureaplasma parvum	●	●	●
Mycoplasma hominis	●	●	●
Yeast-like fungi			
Candida spp.	●	●	●
<i>Candida albicans</i>	●	●	●
Sexually transmitted infections (STI)			
Chlamydia trachomatis	●		●
Mycoplasma genitalium	●		●
Neisseria gonorrhoeae	●		●
Trichomonas vaginalis	●		●
Herpesviruses			
HSV-1	●		●
HSV-2	●		●
CMV	●		●
Human papillomaviruses			
HPV 16	●		
HPV 18	●		
HPV 45	●		
HPV 31/33/35/39/51/52/56/58/59/66/68	●		

How to Choose a Test?

Clinical situation	FB®II	Prima Screen	Secunda Screen
Infectious diseases of reproductive system			
Inflammatory symptoms	●	●	
Symptoms of BV	●		●
Asymptomatic cases			
Preparation for pelvic surgery	●	●	
Preventive screening	●	●	
Treatment monitoring			
Treatment of infectious diseases	●	●	
BV treatment	●		●
Obstetrics			
Pregnancy planning, ART	●	●	
Prophylactic examinations during pregnancy	●	●	
Infertility, miscarriage	●	●	

Biomaterial	Epithelial swab from cervical canal and vagina; cervical swabs in liquid-based cytology transport medium; self-collected vaginal samples
Transport medium	STOR-F, STOR-M; liquid-based cytology media (PreservCyt®, BD SurePath™, EASYPREP, CellPrep)
Tests per kit	12 (24) tests — package S • 24 (48) tests — packages A and A-TL (including control samples)
Storage	+2 °C ... 8 °C (packages S, A, A-TL) • -18 °C ... -22 °C (TechnoTaq MAX polymerase, package A) • 12 months shelf life
Time of analysis	From 2.5 hours (with sample pretreatment)

ANDROBIOME® Line

Androbiome® is a comprehensive real-time PCR-based assay for the assessment of the state of the male urogenital tract microbiota and detection of obligate pathogens — causative agents of sexually transmitted infections (STIs).

Androbiome® Screen — shortened version of the assay, which is intended for screening and differential diagnosis of acute forms of urogenital tract diseases.

Androbiome® — detailed assay for accessing of microbiome structure, which detects a wide range of opportunistic microorganisms. This assay is recommended to prescribe for patients with chronic forms of diseases or complaints of ineffective treatment.

Androbiome® Screen	Androbiome®
Normal microbiota	
<i>Staphylococcus</i> spp.	<i>Staphylococcus</i> spp.
<i>Streptococcus</i> spp.	<i>Streptococcus</i> spp.
<i>Corynebacterium</i> spp.	<i>Corynebacterium</i> spp.
<i>Lactobacillus</i> spp.	<i>Lactobacillus</i> spp.
Opportunistic microorganisms	
<i>Gardnerella vaginalis</i>	<i>Gardnerella vaginalis</i>
<i>Ureaplasma urealyticum</i>	<i>Megasphaera</i> spp. / <i>Veillonella</i> spp. / <i>Dialister</i> spp.
<i>Ureaplasma parvum</i>	<i>Sneathia</i> spp. / <i>Leptotrichia</i> spp. / <i>Fusobacterium</i> spp.
<i>Mycoplasma hominis</i>	<i>Ureaplasma urealyticum</i>
<i>Enterobacteriaceae</i> / <i>Enterococcus</i> spp.	<i>Ureaplasma parvum</i>
<i>Candida</i> spp.	<i>Mycoplasma hominis</i>
	<i>Atopobium</i> cluster
	<i>Bacteroides</i> spp. / <i>Porphyromonas</i> spp. / <i>Prevotella</i> spp.
	<i>Anaerococcus</i> spp.
	<i>Peptostreptococcus</i> spp. / <i>Parvimonas</i> spp.
	<i>Eubacterium</i> spp.
	<i>Haemophilus</i> spp.
	<i>Pseudomonas aeruginosa</i> / <i>Ralstonia</i> spp. / <i>Burkholderia</i> spp.
	<i>Enterobacteriaceae</i> spp. / <i>Enterococcus</i> spp.
	<i>Candida</i> spp.
Obligate pathogens (STIs)	
<i>Mycoplasma genitalium</i>	<i>Mycoplasma genitalium</i>
<i>Trichomonas vaginalis</i>	<i>Trichomonas vaginalis</i>
<i>Neisseria gonorrhoeae</i>	<i>Neisseria gonorrhoeae</i>
<i>Chlamydia trachomatis</i>	<i>Chlamydia trachomatis</i>

ENTEROBIOME® Kids

Reagent kit for the analysis of gut microbiota composition in stool samples from children by real-time PCR

The human gut microbiota is a complex, dynamic ecosystem inhabited by a huge number of different microorganisms, including bacteria, fungi, archaea, and viruses. They interact with each other and with the human organism, participating in a variety of metabolic processes.

The gut microbiota is involved in carbohydrate and protein metabolism, lipid exchange, and also produces important metabolic intermediates such as short-chain fatty acids (SCFAs), secondary bile acids (BAs), vitamins and lipopolysaccharides. These metabolites act as signaling molecules and are involved in the functioning of all systems of the human organism.

Childhood is a period of gut microbiome formation

The establishment of the infant's own gut microbiota is a long and dynamic process in which four successive phases are distinguished. The composition and diversity of the gut microbiota change significantly during childhood, influenced by delivery mode, feeding type, diet, antibiotic use, and gestational age.

Kit format	12 tests per kit (including control samples)
Catalogue No.	M-EBK-0701v01-12
Biomaterial	Stool samples from children
Method	Real-time PCR (multiplex qPCR)
Storage	+2 °C ... 8 °C

HPV Testing Solutions

Real-time PCR diagnostics for human papillomavirus (HPV) detection and genotyping. Three kits offering qualitative screening through comprehensive quantitative analysis with individual genotyping.

QUALITATIVE Genomer HPV Screen 14 HR Types 96 reactions/kit	QUANTITATIVE Genomer HPV Quant 15 15 Types (13 HR + 2 LR) 48 reactions/kit	QUANTITATIVE Genomer HPV Quant 21 21 Types (18 HR + 3 LR) 24 reactions/kit
--	---	---

Key Specifications

Feature	HPV Screen	HPV Quant 15	HPV Quant 21
Method	Qualitative RT-PCR	Quantitative RT-PCR	Quantitative RT-PCR
Reactions/kit	96	48	24
Total genotypes	14	15	21
Sensitivity	5 copies/tube	5 copies/tube	5 copies/tube
Sample types	Cervical/vaginal epithelial swabs	Epithelial swabs, prostate fluid, ejaculate, urine, biopsy	Epithelial swabs, prostate fluid, ejaculate, urine, biopsy
Intended use	Cervical cancer screening	Cervical cancer screening, CIN monitoring	Cervical cancer screening, CIN II/III monitoring
Storage	+2–8 °C	+2–8 °C	+2–8 °C
Shelf life	12 months	12 months	12 months



Genomer HPV Screen

Qualitative Real-Time PCR • 14 High-Risk Types • 96 Reactions/Kit

Qualitative DNA detection of 14 types of human papillomavirus (HPV) with high carcinogenic risk associated with cervical cancer. Provides differential (individual) detection of HPV types 16, 18, and 45; the remaining 11 types are detected as a pooled result.

HPV Types Detected

Individually genotyped	HPV 16, HPV 18, HPV 45
Pooled group detection	HPV 31, 33, 35, 39, 51, 52, 56, 58, 59, 66, 68

Clinical Indications

- Cervical cancer screening and preventive screening programs
- Dynamic monitoring of HPV-associated diseases
- Colposcopic or cytologic abnormalities
- History of other sexually transmitted infections

Technical Specifications

Method	Qualitative real-time PCR (TaqMan hydrolysis probes)
Reactions per kit	96
Tube format	0.2 mL tubes / 8-tube strips
Analytical sensitivity	5 DNA copies per amplification tube
Analysis time	Approx. 1.5 hours (excluding DNA extraction)
Sample types	Epithelial swabs from cervical canal and vaginal mucosa
Detection channels	FAM, HEX, ROX, Cy5, Cy5.5 (5-channel)
Storage	+2 °C to +8 °C • 12 months shelf life

Kit Contents

Ready-to-use PCR mix (preloaded strips) • Internal Control (IC) • Sample Intake Control (SIC) for human DNA assessment • Positive Control (cloned HPV genome) • Negative Control



Genomer HPV Quant 15

Quantitative Real-Time PCR • 15 Types (13 HR + 2 LR) • 48 Reactions/Kit

In vitro DNA test for specific identification and quantification of 15 HPV genotypes. Provides both absolute quantification (log copies/sample) and relative quantification normalized to 10^5 human cells for viral load monitoring.

HPV Types Detected

High-risk group 1	HPV 16, 31, 33, 35, 52, 58 (detected as group)
High-risk group 2	HPV 18, 39, 45, 59 (detected as group)
Individually identified HR	HPV 51, HPV 56, HPV 68
Low-risk (grouped)	HPV 6, HPV 11

Quantification Method

Absolute quantification	Viral load as log (copies/sample) — common logarithm of HPV DNA copy number per sample
Relative quantification	HPV DNA normalized to human cell count (per 10^5 cells) — biomarker for viral replication activity
Sample Intake Control	SIC estimates human DNA amount; SIC < 4.0 indicates insufficient sampling (repeat required)

Technical Specifications

Method	Quantitative real-time PCR (fluorescent TaqMan probes)
Reactions per kit	48
Tube format	0.2 mL × 8-tube PCR strips
Analytical sensitivity	5 copies per amplification tube
Sample types	Epithelial cell swabs, prostate fluid, ejaculate, urine, biopsy
Detection channels	FAM, HEX (IC), ROX, Cy5
Storage	+2 °C to +8 °C • 12 months shelf life

Kit Contents

Preloaded PCR mix (24 × 8-PCR strips) • MAX-Taq polymerase (4 × 500 µL) • Mineral oil (4 × 1.0 mL) • Positive Control (160 µL, cloned HPV genome)



Genomer HPV Quant 21

Quantitative Real-Time PCR • 21 Types (18 HR + 3 LR) • 24 Reactions/Kit

The most comprehensive HPV genotyping kit in the Genomer HPV line. In vitro DNA test for specific identification and quantification of 21 HPV types. Provides absolute and relative quantification with clinically significant viral load thresholds for CIN monitoring.

HPV Types Detected

High-risk types (18)	HPV 16, 18, 26, 31, 33, 35, 39, 45, 51, 52, 53, 56, 58, 59, 66, 68, 73, 82
Low-risk types (3)	HPV 6, 11, 44

Clinical Indications

- Cervical cancer screening (especially women over 30)
- Monitoring therapy for CIN II and CIN III
- ASC-US (atypical squamous cells) confirmatory testing
- HPV diagnostics in oncopathology of cervix, penis, prostate, anus, and larynx
- Inflammatory diseases of the genitourinary system in women and men

Technical Specifications

Method	Quantitative real-time PCR (fluorescent TaqMan probes, hot-start)
Reactions per kit	24
Analytical sensitivity	5 copies per amplification tube (10^3 copies/mL DNA sample)
Diagnostic specificity	99.9% (95% CI: 99.8–99.9%)
Clinically significant load	$>10^3$ copies HPV DNA per 10^5 human cells
Sample types	Epithelial cell swabs, prostate fluid, ejaculate, urine, biopsy
Detection channels	FAM, HEX (IC), ROX, Cy5 (4-channel)
Storage	+2 °C to +8 °C • 12 months shelf life

Kit Contents

Preloaded PCR mix (24 × 8-PCR strips) • MAX-Taq polymerase (4 × 500 µL) • Mineral oil (4 × 1.0 mL) • Positive Control (160 µL, cloned HPV genome, characterized by DNA sequencing)



Advanced Molecular Diagnostics

Biyoteknoloji Arge San. ve Tic. Ltd. Şti.

www.genomer.com.tr

sales@genomer.com.tr